

OAHN Swine Small Scale Herd Postmortem Project

Gastrointestinal Disease in Nursery (Weaner) Pigs (<20 kg, 1-3 animals)

TISSUE SAMPLING LIST - from EACH pig:

Fresh tissues: 1 sample per Whirl-Pak bag, per pig – label with stickers provided.

Formalin-fixed tissues: 1 sample jar per pig.

Tissue	Fresh (<i>include lesions</i>)	In Formalin (<i>include lesions</i>)
Tonsil	X	X
Thyroid gland		X
Trachea		X
Thymus		X
Lung	X 2	X
Heart		X
Liver	X	X
Spleen	X	X
Kidney	X	X
Adrenal gland		X
Pancreas		X
Mesenteric lymph node	X	X
Gastrohepatic lymph node	X	X
Stomach		X
Duodenum		X
Jejunum		X
Ileum	X 4	X
Spiral colon	X 3	X
Feces	X	
Diaphragm		X
Skeletal muscle		X
Brain (or head)	X half, sagittally sectioned (OR entire head, fresh)	X half, sagittally sectioned (OR entire head, fresh)



**ANIMAL HEALTH
LABORATORY**

Gastrointestinal Disease in Nursery (Weaner) Pigs (<20 kg, 1-3 animals)

Fresh tissues: Use a separate Whirl-Pak bag for each sample – use provided stickers and marker to label bags with the tissue name and pig designation (pig A, B, or C).

Formalin-fixed tissues: Use 1 sample jar per pig and label with pig designation (pig A, B, or C). For solid organs, samples should be approximately 1-2 cm². For intestine, take approximately 3 cm long segments and partially open longitudinally with scissors (to allow formalin to contact the mucosal surface).

Test list (for AHL use - included on submission form):

Virology: PRRSV PCR *prtqn*, influenza A virus PCR *inflpqr*: lung (pool)

Rotavirus A/B/C PCR *rotapqr*, porcine coronavirus PCR *pcorpcr*: ileum (pool)

Bacterial culture *bsetup*, *cults*: ileum, colon (pool by pig)

Lawsonia PCR *lapqr*: ileum (pool)

Brachyspira spp PCR *brpqr* and *bspqr*: colon (pool)

Parasitology – fecal flotation *fflot*: feces (pool)

Histopathology *hist*

Freeze and hold: tonsil / spleen, lung, kidney / liver, lymph nodes, ileum, colon, brain.



ANIMAL HEALTH
LABORATORY