

**OAHN Swine Small Scale Herd Postmortem Project**

**Gastrointestinal Disease in Sows /Gilts / Boars**

**TISSUE SAMPLING LIST**

**Fresh tissues:** 1 sample per Whirl-Pak bag, per pig – label with stickers provided.

**Formalin-fixed tissues:** 1 sample jar per pig

<b>Tissue</b>	<b>Fresh (<i>include lesions</i>)</b>	<b>In Formalin (<i>include lesions</i>)</b>
Tonsil	X	X
Thyroid gland		X
Trachea		X
Thymus		X
Lung	X 2	X
Heart		X
Liver	X	X
Spleen	X	X
Kidney	X	X
Adrenal gland		X
Pancreas		X
Mesenteric lymph node	X	X
Gastrohepatic lymph node	X	X
Stomach		X
Duodenum		X
Jejunum		X
Ileum	X 4	X
Spiral colon	X 3	X
Feces	X	
Diaphragm		X
Skeletal muscle		X



**ANIMAL HEALTH  
LABORATORY**

Brain (or head)	X half, sagittally sectioned (OR entire head, fresh)	X half sagittally sectioned (OR entire head, fresh)
--------------------	---	--

## Gastrointestinal Disease in Sows /Gilts / Boars

**Fresh tissues:** Use a separate Whirl-Pak bag for each sample – label bags with the stickers provided (tissue name).

**Formalin-fixed tissues:** Use 1 sample jar per pig. For solid organs, samples should be approximately 1-2 cm<sup>2</sup>. For intestine, take approximately 3 cm long segments and partially open longitudinally with scissors (to allow formalin to contact the mucosal surface).

### Test list (for AHL use - included on submission form):

Virology: PRRSV PCR *prrsv*, influenza A virus PCR *inflpcr*: lung

Rotavirus A/B/C and Porcine coronavirus PCR *rotapcr*, *pcorpcr*: ileum

Bacterial culture *bsetup*, *cults*: ileum, colon (pool)

Lawsonia PCR *lapcr*: ileum

Brachyspira spp PCR *brpcr* and *bspccr*: colon

Parasitology – fecal flotation *fflot*: feces

Histopathology *hist*

Freeze and hold: tonsil / spleen, lung, kidney / liver, lymph nodes, ileum, colon, brain.



ANIMAL HEALTH  
LABORATORY