

OAHN Swine Small Scale Herd Postmortem Project

Lameness in Nursery (Weaner) Pigs (1-3 animals <20 kg)

TISSUE SAMPLING LIST - Collect the following samples from EACH pig:

Fresh tissues: 1 sample per Whirl-Pak bag, per pig – label with stickers provided.

Formalin-fixed tissues: 1 sample jar per pig

Tissue	Fresh (include lesions)	In Formalin (include lesions)
Tonsil	X	X
Thyroid gland		X
Trachea		X
Thymus		X
Lung	X 2	X
Heart		X
Liver	X	X
Spleen	X	X
Kidney	X	X
Adrenal gland		X
Pancreas		X
Mesenteric lymph node	X	X
Gastrohepatic lymph node	X	X
Stomach		X
Duodenum		X
Jejunum		X
Ileum		X
Spiral colon		X
Diaphragm		X
Skeletal muscle		X
Joints: 1 representative joint with lesions, per pig	X synovium / joint capsule, swab, or aspirate X2	X synovium / joint capsule
Brain (OR submit head intact)	X half, sagittally sectioned	X half, sagittally sectioned
Spinal cord (OR vertebral column)	X	X



**ANIMAL HEALTH
LABORATORY**

Lameness in Nursery (Weaner) Pigs (1-3 animals <20 kg)

Fresh tissues: Use a separate Whirl-Pak bag for each sample – use stickers and marker provided to label bags with the tissue name and pig designation (pig A, B, or C).

Formalin-fixed tissues: Use 1 sample jar per pig and label with pig designation (pig A, B, or C). For solid organs, samples should be approximately 1-2 cm². For intestine, take approximately 3 cm long segments and partially open longitudinally with scissors (to allow formalin to contact the mucosal surface).

Test list (for AHL use - included on submission form):

Virology: PRRSV PCR *prtqn*, influenza A virus PCR *inflpcr*: lung (pool)

Bacterial culture *bsetup*, *cults*: joint swab(s), joint capsule, or aspirate (pool by pig)

Mycoplasmology: M.hyosynoviae and M.hyorhinis PCR *hsypcr*, *hrhpcr*: joint swab(s), joint capsule, or aspirate (pool by pig)

Histopathology *hist*

Freeze and hold: tonsil / spleen, lung, kidney / liver, lymph nodes, brain, spinal cord.



ANIMAL HEALTH
LABORATORY