

OAHN Swine Small Scale Herd Postmortem Project

Respiratory Disease in Nursery (Weaner) Pigs (<20 kg, 1-3 animals)

TISSUE SAMPLING LIST

Collect the following samples from EACH pig:

Fresh tissues: 1 sample per Whirl-Pak bag, per pig – label with stickers provided.
Formalin-fixed tissues: 1 sample jar per pig.

Tissue	Fresh (<i>include lesions</i>)	In Formalin (<i>include lesions</i>)
Tonsil	X	X
Thyroid gland		X
Trachea		X
Thymus		X
Lung	X 4	X
Heart		X
Liver	X	X
Spleen	X	X
Kidney	X	X
Adrenal gland		X
Pancreas		X
Mesenteric lymph node	X	X
Gastrohepatic lymph node	X	X
Stomach		X
Duodenum		X
Jejunum		X
Ileum		X
Spiral colon		X
Diaphragm		X
Skeletal muscle		X
Brain (or head)	X half, sagittally sectioned (OR entire head, fresh)	X half, sagittally sectioned (OR entire head, fresh)



**ANIMAL HEALTH
LABORATORY**

Respiratory Disease in Nursery (Weaner) Pigs (<20 kg, 1-3 animals)

Fresh tissues: Use a separate Whirl-Pak bag for each sample – use provided stickers and marker to label the bag with the tissue name and pig designation (pig A, B, or C).

Formalin-fixed tissues: Use 1 sample jar per pig and label with pig designation (pig A, B, or C). For solid organs, samples should be approximately 1-2 cm². For intestine, take approximately 3 cm long segments and partially open longitudinally with scissors (to allow formalin to contact the mucosal surface).

Test list (for AHL use - included on submission form):

Virology: PRRSV PCR *prtqn*, influenza A virus PCR *inflpcr*: lung (pool)

Mycoplasmology: *M.hypopneumoniae* PCR *mhpcr*: lung (pool)

Bacterial culture *bsetup, cults*: lung (pool).

Histopathology *hist*

Freeze and hold: tonsil /spleen, lung, kidney /liver, lymph nodes, brain.



ANIMAL HEALTH
LABORATORY